## NL

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## ENG

**Indications:**

Control of knee mobility. Post-surgical knee treatment. Rehabilitation after injury.

**Precautions:**

- Only a doctor is qualified to prescribe the orthosis and determine the duration of treatment.
- Read these instructions carefully before using the orthosis. If you have any questions, consult your doctor or the shop where you purchased it.
- The orthosis requires assembly and prior adjustment by qualified personnel.
- Although the orthosis is not for single use, it should only be used by one patient. The product should not be reused by other patients.
- The general instructions of use and specific indications determined by your doctor should always be followed.
- Treatment with the orthosis should be regularly monitored by the doctor.
- If any side effects are experienced, your doctor should immediately be informed.
- Your doctor should be informed if a skin condition appears or any type of sensitisation occurs.
- Your doctor might consider the use of topical creams in conjunction with the orthosis.
- In areas of support where pressure is applied, the skin should not be damaged or hypersensitive.
- For the orthosis to function properly, its components must be regularly checked to ensure that they are in good condition.
- Do not let children play with this device.
- Make sure that the Velcro fasteners are free from fluff and lint to ensure that their closure properties are preserved.
- Clean the product regularly and maintain proper hygiene.
- The medical staff responsible for monitoring treatment may explain the suitability of the product to the patient or suggest its replacement if it is damaged or worn.
- Do not wear the orthosis while having a bath or shower, or sleeping at night.

**Instructions for proper fitting:**

These knee braces are equipped with ROM joints which limit flexion and extension movements by means of interchangeable wedges. They feature fastening straps marked with a blue dot for the proximal area (thigh) and a red dot for the distal area. All models are available in 5 sizes (S, M, L, XL, XXL).

Fig.1: When inserting the leg, the thigh area should be open to make fitting easier. Open the velcro fasteners and fit the brace to the leg. Position the knee brace ensuring that the opening is over the kneecap (A). The joint should enable flexion and be centred on the axis of the knee (B).

Fig.2: Then close the top part of the knee brace (C).

Fig.3: And close the top strap marked with the blue dot (D).

Fig.4: Repeat the operation with the bottom strap marked with the red dot (E). The joint should be adjusted by bending the bars and always using the right tool (Allen key). Adjustment should be made according to the patient's anatomy. Do not handle the inside of the joint where the hinge is located.

The ROM joint can be adjusted by increments of 10°. Range 0° to 80°. There are eight sets of stops interchangeable by means of a screw for both extension and flexion.

**Product maintenance - washing instructions:**

This product is latex free. Wash at 40°C. Do not dry clean. Do not iron. Do not use bleach. Do not tumble dry.

Remove the ROM joints before washing. Do not change the limit stops set by your doctor or orthopaedic surgeon. Clean with a damp cloth. Do not grease the hinges.

**Classification:**

Class I medical device without a measuring function and supplied in a non-sterile condition. This orthosis complies with Council Directive 93/42/EEC, amended by 2007/47/EC, concerning medical devices.